









Ion Chamber Compatibility Chart					
Dosimeters Ion Chambers	2026	9010 / 9015	Accu-Dose 2086 / 2186	Accu-Pro 9096	Accu-Gold AGDM
10X5 Series		✓	10A96* 	10A96* 	10A96* 
10X6 Series			✓	✓	✓
10X9 Series			95A96** 	95A96** 	95A96** 
20X6 Series	✓				

✓ = Direct Connect

* = Adapter Model 10A96 Required 

** = Adapter Model 95A96 Required 

Note: Discontinued model 1015/1515 used 10X5 chambers and discontinued model 9095 used 10X9 chambers.